

Breezy Style X



Status: 11_2021 to 12_2022
Rev:09_2021

ORDER FORM

Errors and imissions excepted.

Billing Address Or Customer-ID			Shipping Address	Identical To Billing Address
C.-ID:			C.-ID:	
Name:			Name:	
Street:			Street:	
Town/PC:			Town/PC:	
E-Mail:			E-Mail:	
Phone:			Phone:	
Order Information			User Weight	
Order	Calculation	Built-4-Me [®]	User Weight Capacity: 125 kg	
Purchase Order No.:			User Weight: kg	
Art. No.: 08000101		Reimbursement Code: xxx		
Order Date:	Day	Month	Year	

Base Price **€ X.XXX,-**

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat Width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10-20 mm altogether).
- b) On the wheelchair: From outside frame tube to outside frame tube.

Seat Depth (SD):

- a) For the user: When sitting upright: leave space for 2-3 fingers between front edge of the seat upholstery and back of knee (approx. 30-60 mm). Care must taken with adj. upholstery.
- b) On the wheelchair: From back tube to front seat upholstery.

Front Seat Height:

- a) For the user: When sitting with 90°knee angle: from back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground. Min. 30 mm longer than user measurement. Remember to include the cushion height.

Rear Seat Height:

On the Wheelchair: From back of seat tube to the ground.

Note: Max. seat height difference rear to front is 110 mm

Back Height:

- a) For the user: From the seat, including cushion, to be no more than 10-20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.

Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably

Back angle:

Important: The back height depends on the degree of disability. 0° is the angle 90° to the seat. Negative angles are <90° and positive angles are >90°.

Centre of Gravity:

On the wheelchair: From the back tube to the middle of the axle plate.

Technical Data

Max. Overall Height:	995 mm
Max. Overall Length:	1140 mm
Overall Width:	SW + 180 mm
Weight	from 14.7 kg

Product is crash tested according to ISO 7176-19 and therefore approved for transportation in a vehicle in combination with crash test option. For transport without a user, fix the chair safely in the vehicle (e.g. car or plane).

Seat Width						
	SW38	SW405	SW43	SW455	SW48	SW52
mm	380	405	430	455	480	520

Seat Depth		€
SD435	435 mm	-

Front Seat Height		€
FSH485	485 mm (standard for 6" and 8" castor wheels)	-

Seat		€
Seat Upholstery		
BRU020065	Seat sling, Nylon, padded, black	-
Cushion		
BRU020066	Seat cushion, foam (medium soft), cover black, 50 mm	XXX,-

Backrest		€
BRU030122	Backrest, fixed, push handles long, lumbar curve 10° (1)	-
BRU030014	Backrest, recline (0° - 30°), push handles long (2)	XXX,-
Backrest Accessories		
BRU030080	Headrest, Skay, height- and angle-adjustable, 260 x 130 mm	XXX,-
Backrest Upholstery		
BRU030076	Backrest sling, Nylon, adjustable, 3 straps, black	-

Armrest		€
BRU040052	Single post sideguard, armpad (330 mm), tool height-, depth-adjustable, removable	-

Legrest & Footplates		€
Hanger Angle		
BRU050002	Hanger angle 70°	-
Composite Hanger		
BRU050003	Hanger, composite, swing away, left	-
BRU050004	Hanger, composite, swing away, right	-
Elevating Legrest (ELR)		
BRU050009	Hanger, aluminium, flip-up (0° - 90°), swing away, left	XXX,-
BRU050010	Hanger, aluminium, flip-up (0° - 90°), swing away, right	XXX,-
Amputee Support Pad		
BRU050014	Amputee support, left (3)	XXX,-
BRU050013	Amputee support, right (3)	XXX,-
Footplate		
BRU050104	Footboard, divided, composite, flip-up (sideways), heel loop, black	-

Castor Wheels		€
BRU080206	Castor, 6", solid, rim: composite, black	-
BRU080049	Castor, 8" x 1 1/4", solid, rim: composite, black	-
BRU080047	Castor, 8" x 2", solid, rim: composite, black	-

Rear Wheels		€
BRU070050	Rear wheel, 24", solid, 36 spokes (silver), cross-wise spokes, handrim Alu anodized-silver	-

(1) Back height 425 mm, adjustable to 400 and 450 mm, COG 25 mm
(2) Standard: BH 480 mm, Backrest sling, Nylon, 4 Straps, height-adj. armrest, stabilizer bar, anti-tip tubes required, COG +25 mm
(3) Not with footplate

Backrest & Accessories



// BRU030122 Backrest, fixed, push handles long, lumbar curve 10°



// BRU030014 Backrest, recline (0° - 30°), push handles long



// BRU030080 Headrest, Sky, height- and angle-adjustable, 260 x 130 mm

Armrest



// BRU040052 Single post sideguard, armpad (330 mm), tool height-, depth-adjustable, removable

Legrest & Footplates



// BRU050002 Hanger angle 70 °



// BRU050009 Hanger, aluminium, flip-up (0° - 90°), swing away, left (ELR)



// BRU050014 Amputee support, left

Rear Wheels



// BRU070050 Rear wheel, 24", solid, 36 spokes (silver), cross-wise spokes, handrim Alu anodized-silver

Axles		€
BRU070150	Quick release axle	-
Center Of Gravity (COG)		€
COG25	Centre of gravity +25 mm	-
Wheel Locks		€
BRU060002	Knee-lever brake	-
BRU060010	Brake lever extension, knee-lever brake, fold-down, pair	XXX,-
Trays		€
BRU090011	Tray, push-fit (1)	XXX,-
Accessories		€
BRU090039	Crutch holder, left (2)	XXX,-
BRU090060	Crutch holder, right (2)	XXX,-
BRU070831	Spoke guard, transparent	XXX,-
BRU090044	Anti-tip tubes, plug-in, flip-up, pair (3)	XXX,-
BRU090042	Anti-tip tubes, plug-in, flip-up, left (3)	XXX,-
BRU090043	Anti-tip tubes, plug-in, flip-up, right (3)	XXX,-
BRU090135	Positioning belt, carseat buckle	XXX,-
BRU090102	Crash test option	XXX,-
Colours		€
BRU100140	Apple green	-
BRU100116	Brilliant red	-
BRU100141	Navy Blue	-
BRU100142	Selenium grey	-
BRU100124	Signal white	-

(1) SW 380 mm to 480 mm

(2) Not with reclining backrest

(3) Standard with reclining backrest

Wheel Locks



// BRU060002 Knee-lever brake



// BRU060010 Brake lever extension, knee-lever brake, fold-down

Trays & Accessories



// BRU090011 Tray, push-fit



// BRU090039/60 Crutch holder, left/right



// BRU070831 Spoke guard, transparent



// BRU090042-43 Anti-tip tubes, plug-in, flip-up, left/right



// BRU090135 Positioning belt, car-seat buckle

Comments for Customisations



Breezy Style X

€

Base price	<input type="text"/>
Total Of Additional Equipment	<input type="text"/>
Packaging And Transport	<input type="text"/>
Exclusive VAT	<input type="text"/>
VAT	<input type="text"/>
Total price	<input type="text"/>

Date

Signature

Day Month Year

**Sunrise Medical S.L**

Poligono Bakiola, 41
 48494 Arrankudiaga
 Phone: +34 902142 434
 Fax: +34 946 481 575
 E-Mail: Monica.garcia@sunrisemedical.es
 Patricia.delafuente@sunrisemedical.es
 www.SunriseMedical.eu