





KRYPTON R

Edition: 11.21 Valid until: 12.22 Rev. Sept. 21

Deliver to: Name: Street:	Deliver to: Name: Street:	
Town: Postal code:	Town:	Postal code:
E-Mail:	E-Mail:	
Tel: Fax:	Tel:	Fax:
Order date: Purchase order: Number required:	Marked for:);
•	• = Standard ne (longer leadtime)	□ = Cost option o = No cost option
Art. no.: 770 00 001	Max. user wei USER WEIGH	·

Base Price

Mandatory information to ensure correct configuration

KRYPTON R: STANDARD CONFIGURATION

KRYPTON R

- · Rigid wheelchair
- · Carbon fibre frame
- · Adjustable back angle
- · Foldable backrest
- · Adjustable COG
- Rear seat height adjustable +/- 1cm
- · Carbon fibre axle-stems
- · Carbon fibre axle tube





Krypton R

= Standard

O = No cost option

□ = Cost option

CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH SELECTION - in mm; measured from outside frame tube to outside frame tube

320	340	360	380	400	420	440	460
0	0	0	0	0	0	0	0

SEAT DEPTH SELECTION - in mm

340	360	380	400	420	440	460	480
0	0	0	0	0	0	0	0

FRONT SEAT HEIGHT SELECTION - in dependance of front castor size, fork, frame height version The high frame is 40 mm higher than the low frame.

Front seat height in mm

							_	_	_		_		_	_		
_		Fork (length)	Frame:	430	440	450	460	470	480	490	500	510	520	530	540	550
		111 mm		0	0	0										
		123 mm	low			0	0	0								
		Frog Leg	low													
	3"	One Arm														
		111 mm						0	0	0						
		123 mm	high							0	0	0				
		Frog Leg	nign													
		One Arm														
		111 mm			0	0	0	0								
		123 mm														
	4"	Frog Leg	low													
		One Arm														
ē		111 mm							0	0	0	0				
ğ		123 mm	himb													
Ę		Frog Leg	high													
Front Caster		One Arm														
- Г		111 mm					0	0	0							
		123 mm	low													
		Frog Leg	low													
		One Arm														
	5"	111 mm									0	0	0			
		123 mm	himb													
		Frog Leg	high													
		One Arm	1													
	e"	123 mm	himb													
	6"	Frog Leg	high													

REAR SEAT HEIGHT SELECTION - in dependance of rear wheel size

The rear seat height is adjustable ± 10 mm.

The rear seat height can only be between 0 cm and 110 mm below the selected front seat height.

Caution: with rear seat heights 460 mm and above, the width between the sideguards increases by 35 mm

Rear seat height in mm

No upwards-adjustment of rear-seat-height possible. Only -1 cm adjustable.

Rear Wheels	380	390	400	410	420	430	440	450	460	470	480	490	500
24"	O ⁽¹⁾	$O^{(2)}$	$O^{(2)}$	0	0	0	0	0	0	0	0	$O^{(3)}$	
25"		O ⁽¹⁾	O ⁽²⁾	O ⁽²⁾	0	0	0	0	0	0	0	0	O ⁽³⁾

(1):	No adjustment of rear seat-height possible
(2):	No downwards-adjustment of rear-seat-height possible. Only +1 cm adjustable.

For more information www.SunriseMedical.eu © Sunrise Medical 2021 & 22

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

a) For the user: Allow space between sideguard and hip on each side of the chair

(approx. 10 - 20 mm altogether)

b) On the wheelchair: From outside frame tube to outside frame tube.

Seat depth (SD):

a) For the user: When sitting upright: leave space for 2 - 3 fingers between

front edge of the seat upholstery and back of knee

(approx. 30 - 60 mm). Care must be taken with adj. upholstery.

b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

a) For the user: When sitting with 90° kneeangle:

From back of knee to bottom of shoe sole.

b) On the wheelchair: From front of seat tube to the ground .

(min. 30 mm longer than user measurement) Remember to include the cushion height.

Rear seat height:

b) On the wheelchair: From back of seat tube to the ground.

Note: Max. seat height difference rear to front is 130 mm

Back height:

a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.

b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.

Note: The user's adility to sit upright has a major influence on the back height

and it can therefore differ considerably


Back angle: Important: The back height depends on the degree of disability.

 0° is the angle 90° to the seat. Negative angles are ${<}90^{\circ}$

and positive angles are >90°.

Low Seat Frame vs. High Seat Frame:

The high frame is 40 mm higher than the low frame. So you can choose shorter and lighter forks to minimize the weight of your wheelchair.



Technical Data Overview:

Max. overall height: 950 mm Max. overall width: 24" rear wheel 25" rear wheel Camber Max. overall length: 930 mm For 0° camber: SW + 170 mm SW + 170 mm Max. turning circle: 700 mm For 3°camber: SW + 230 mm SW + 240 mm Max. safe slope: 7°

Wheelchair weight: from 6,2 kg
Dismantled weight: from 3,7 kg
Weight of heaviest part: 2,1 kg

All changes in the configuration might change the seat heights.

Manufacturing tolerances are +/- 10 mm for all measurements.

© Sunrise Medical 2021 & 22

WEIGHT

We want to be very transparent on weight. The following information will be a guide to the lightest possible Krypton R: how to achieve a 6,2 kg wheelchair, and how to order the lightest possible configuration for your Krypton R

Krypton R: 6,2 kg configuration Transport weight: 3,7 kg

Frame: high frame, 92°, 30 mm inset, seat width: 360 mm, seat depth: 400 mm. Seat height: front 470 mm, rear 420 mm. Seat-sling Airflow superlight.

Wheel-lock: lightweight compact wheel-lock

Footrest: tubular footrest



Backrest: Carbon backrest, EXO superlight backrest-upholstery, without push-handles

Sideguards: Carbon fibre, slim-style

Rear wheels: 24" Proton wheels with Aluminium handrims (coated or anodized) and Schwalbe One tyres

Forks and casters: 111 mm carbotecture fork, 3" caster

Please mind that there is a 5% tolerance on all weight specifications due to weight deviations on mounting parts (tyres, caster-wheels etc.)

LIGHTWEIGHT-GUIDE FOR ORDERING THE KRYPTON R

Sunrise Medical offers a wide variety of options for the Krypton R to suit every user's needs. However, most options add significant weight to the wheelchair. To make this more transparent when filling in the order-form, and to assist you in ordering the lightest possible wheelchair, we indicate the approximate weight-increase of the options on the order form:

- **(L)** indicates the lightest possible option within the range of options
- (+) indicates that this option adds more than 500 g to the wheelchair's weight compared to the lightest option within the range

LIGHTWEIGHT OPTIONS FOR THE KRYPTON R

We developed several options for the Krypton R to save weight on the chair without sacrifing daily-use suitability, comfort or durability



Carbon backrest (KRR030006): full carbon backrest with very close cross-tube (20mm crosstube to backrest-upholstery) (-100 g compared to Standard backrest)



Seat sling "airflow" superlight (KRR020019): fully adjustable seat-sling with composite buckles and lighter straps (max. user weight 85 kg) (-100 g compared to seat sling "airflow")



EXO superlight (KRR030005): lightweight version of the EXO backrest upholstery. EVA-foam padding, adjustable straps, wider cut of the upholstery to cover backrest tubes.

(-160 g compared to Std. EXO backrest-upholstery)



Axle-plug superlight (KRR070049): (only 0° camber): glued into the axle-tube, no adjustability of the track-width possible (-160 g compared to widthadjustable

axle-plug KRR070600 / 603)

Krypton	R		
	•	= Standard O = No cost option \square = Cost option	Price in €
		SEAT UPHOLSTERIES AND CUSHIONS	
Seat Sling			
KRR020012	0	Seat sling "airflow" with 1 utility bag	
KRR020017		Seat sling "airflow" with 2 utility bags Seat sling "airflow" with utility bag and catheter bag (from SW 340)	21
KRR020018 KRR020019		Seat sling "airflow" with utility bag and catheter bag (from SW 340 Seat sling "airflow" superlight. Without bags, max user weight 85kg (1) (L	
KRR020019 KRR020000	0	Without seat sling	-)
Seat cushion Foam: (mediu			
Latex: (very s		ft) Thickness: KRR020101 □ 30 mm KRR020102 □ 50 mm Thickness: KRR020103 □ 50 mm	
		or complete product range see Jay order forms	
JAY020002 JAY020003		JAY Soft Combi P	-
JAY020003 JAY020001		JAY Easy Visco, Curved Base	- -
JAY020007		JAY Easy Fluid, Curved Base	
JAY020004		JAY Lite	
_	-=:00		
Seat accesso KRR090018	ories	Positioning belt (only with Standard backrest)	
KKKU9UU10		Positioning belt (only with standard backlest)	
		BACKREST	
KRR030001	0	Standard backrest (aluminium, black)	
		Choose back height factory setting (mm). 250 - 300 300 - 350 350 - 400 400 - 450	
		Height adjustments possible in the O 250 O 300 O 350 O 400	
		range 250-300; 300-400; 400-450 mm. O 275 O 325 O 375 O 425	
		O 300 O 350 O 400 O 450	
KDD000006		Only the least fined being the constitute of the unbelaton. "Vanton" mikhor grin (2) (1)	
KRR030006		Carbon backrest - fixed height, crosstube very close to the upholstery, "Krypton" rubber grip (2) (L) Choose back height (mm)	
		no height adjustment possible. O 250 O 325 O 400	
		O 275 O 350 O 425	
		O 300 O 375 O 450	
		0 300 0 370 0 400	
KRR030012	•	Backrest, angle-adjustable (non foldable) (-31° to +15°) (only in black anodizing)	
KRR030015	0	Backrest, angle-adjustable, fold-down (-31° to +15°) (only in black anodizing)	
KRR030102		Backrest, angle-adjustable, fold-down, locked when folded, colour: black (-31° to +15°)	
		colour: □ Orange □ Silver □ Red □ Blue	
		KRR030036 KRR030037 KRR030038 KRR030039	
Pushhandles	s: chr	pose only one option	
KRR030200	0	Omit push handles (L	_)
KRR030201	0	Push handles long	
KRR030203		Fold-down push handles (from 325mm backrest-he	eight)
KRR030204		Height adj. push handles (+) (only with Standard backrest; from 325mm backrest he	-
KRR030205		Height adj. push handles for Jay backs (+) (only with Standard backrest; from 325mm backrest he	eigth)
Backrest uph	10 st	ery	
KRR030316	0	EXO Backrest sling, padded, adjustable, 1 accessory bag, black (EXO Evo version)	
		Colour, binding tape: ● Black O Orange O Silver O Red O Blue	
KRR030317		EXO PRO Backrest, padded, adj., breathable spacer fabric, 1 accessory bag, black (EXO Evo version	4)
		Colour, binding tape: Black O Orange O Silver O Red O Blue	
KRR030005		EXO superlight, padded, adjustable, black (EXO Evo version) (L	-)
**************************************		Colour, binding tape: Black O Orange O Silver O Red O Blue Without backrost unbelstory	
KRR030300	0	Without backrest upholstery	
(1):	Мах.	user weight 85 kg, only available up to seat-width 400 mm	For more information
(2):	Not	with height adj push handles, not with Jay J3 backs, not with positioning-belt	www.SunriseMedical.eu
	(L)=	lightest option (+)= adds 500 g or more to the wheelchair-weight	Sunrise Medical 2021 & 22

SEAT UPHOLSTERIES AND CUSHIONS



KRR020012 Seat sling, airflow



KRR020019 Seat sling, airflow, superlight

BACKREST



KRR030006 Carbon Backrest



KRR030102 Backrest, angle-adjustable, lock when folded



KRR030201 Push handles long



KRR030203 Push handles, fold-down



KRR030204 Push handles , height-adjustable



KRF030316 EXO backrest sling (EXO Evo



KRF030317 EXO PRO backrest sling (EXO Evo version)



KRF030005 EXO superlight backrest sling (EXO Evo version)

FOOTREST AND FOOTPLATES



KRR050029 Tubular footrest lightweight



KRR050025 Platform footrest LW composite plate



KRR050026 Platform footrest LW carbon plate



KRR050125 Performance footrest alu angle adj



KRR050127 Foot positioning plate

Krypton I	R		
	•	·	rice in €
		JAY 3 BACKS AND ACCESSORIES - only with Standard backrest (+)	
Jay 3 Back -			m
Mid Thoracic: Upper Thorac		□ JAY030012 Short - 420 mm □ JAY030019 Medium - 340 mm □ JAY030033 Tall - 380 m □ JAY030012 Short - 420 mm □ JAY030019 Medium - 460 mm □ JAY030025 Tall - 500 m	
Shoulder Heigh		□ JAY030015 Short - 530 mm □ JAY030022 Medium - 570 mm □ JAY030028 Tall - 610 m	
Jay 3 Back - Mid Thoracic:		□ JAY030011 Medium - 340 mm □ JAY030018 Tall - 380 m	m
Upper Thorac		□ JAY030013 Short - 420 mm □ JAY030020 Medium - 460 mm □ JAY030026 Tall - 500 m	
Shoulder Heigh		□ JAY030016 Short - 530 mm □ JAY030023 Medium - 570 mm □ JAY030029 Tall - 610 m	
		k - Shallow Contour	
Lower Thorac		□ JAY060001 Short - 170 mm □ JAY060002 Tall - 240 m	m
Mid Thoratic:	,io.	□ JAY060003 Short - 300 mm □ JAY060004 Tall - 380 m	
	4 1	AV 12 Peaks	
Accessories WHT030415	tor J		3 hacks)
WHT030416		□ Whitmyer Headrest, COBRA XTRA, flip back (only with Jay 3 □ Whitmyer Headrest, AXYS, Standard (only with Jay 3	
JAY090058		☐ JAY Accessory bag (only with Jay 3 Carbon backs)	, 5401.0)
		(c)	
		SIDEGUARDS AND ARMRESTS	
KRR040101			L)
KRR040102		Sideguard, aluminium, cold weather protection (colour: black, matt)	
KRR040103		Sideguard, aluminium, fender, cold weather protection (colour: black, matt) (min rear seat height 3	
KRR040104 KRR040107		Sideguard, carbon fibre, fender, cold weather protection, black (min rear seat height 3	(90 mm)
KRR040052		Sideguard, composite, removable, 21 cm high (31 cm long), black Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable (2) (+)	
KRR040117		Tubular armrest, Alu, swing away (sideways), padded, removable (+)	
KRR040000	0	w/o Sideguard	
		FOOTREST AND FOOTPLATES	
KRR050029	•	, , , , , , , , , , , , , , , , , , , ,	L)
KRR050030		Footrest, tubular, lightweight, aluminium, with cover	
KRR050025 KRR050024		Platform footplate lightweight angle & depth adjustable composite	
KRR050024		Platform footplate lightweight angle & depth adjustable aluminium (colour: black, matt) Platform footplate lightweight angle & depth adjustable carbon	
KRR050125		Performance footplate angle adjustable (3)	
KRR050128		Footboard, platform, plastic, angle and depth adjustable, flip-up (to the back) (4)	
KRR050127		Adjustable side protection, to create a side-stop on the footboard (5)	
		LOWER LEG LENGTH (LLL) Attention: LLL has to be 30 mm less than seat height front; all mea	asurements in mm
	0	360 O 370 O 380 O 390 O 400 O 410 O 420 O 430	
	0	440 O 450 O 460 O 470 O 480 O 490 O 500 O 510	
		FRAME VERSIONS	
1/DD010000		Frame length	
KRR010206 KRR010207	0	Frame standard Frame long (20 mm additional frame length) (not in combination with SD 480 mm)	
KKK010207			
		Frame angle	
KRR010012	•	100° frame	
KRR010013		92° Performance frame	
		Frame inset	
KRR010101	•	Inset: Small (10 mm)	
KRR010103		Inset: Performance (30 mm) (not available for seat width 320 mm)	
		Castor position (see explanation on next page)	
KRR010301	•	Castor house position wide	
KRR010302	0	Castor house position narrow (not with seat-depths < 400mm, please see recommendations on next pa	age)
(4)	1454	2// 1 1 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	
(1) (2)		24" rear wheel min. rear seat height 410 mm; with 25" rear wheel min. rear seat height 430 mm with COG 56 mm	For more info
(2)	_	with COG 56 mm	For more information www.SunriseMedical.eu
(3) (4)		weight limitation 75 kg, not with 4" castor	© Sunrise Medical 2021 & 22
(4) (5)		for platform footplates lightweight angle & depth adjustable	C Carriot modical Evel & ZZ
1.7	_	ghtest option (+)= adds 500 g or more to the wheelchair-weight	

SIDEGUARD & ARMREST



KRR040101 Slim style



KRR040104 Carbon fibre, cold weather protection, fender



KRR040107 Sideguard, composite, removable (black)




KRR040117 Tubular armrest, swing away, removable, padded




KRR040052 Single-post sideguard, armpad (33 cm)

FRAME VERSIONS: CASTOR HOUSE POSITIONS



Further to the front and wider position of the casters for optimized rolling resistance

Position narrow:





22 mm further rearwards and 24 mm per side more narrow position of the casters for a more narrow front-end and improved manoeuverability in tight spaces. Reduced front-stability compared to <u>position wide</u>, we recommend this option in combination with seat-height differences (front to rear) of 40mm or more or in combination with KRR010207 "frame, long".



KRR010301 Castor house position wide



KRR010302 Castor house position narrow



Castor house position wide vs. narrow

Kryptor	ı R		
	•	= Standard O = No cost option	Price in €
		CENTRE OF GRAVITY (COG) - all measurements are in mm; measured from front of back post; choose only	y one option
COG56	0	56 mm ° 0	
COG68	0	68 mm	
COG80	0	80 mm	
COG92	0	92 mm	
COG104	0	104 mm	
		FRAME ACCESSORIES	
KRR090031		Suspension up to 70 kg (1) (+) (chair will be built using aluminium axle-stems and axle-tube)	
KRR090032		Suspension up to 120 kg (1) (+) (chair will be built using aluminium axle-stems and axle-tube)	
		Tip assist (step tube) (2)	
KRR090001		Left KRR090002 □ Right	
KRR090003		Pair (+) (not available with anti-tip)	
KRR090009		Crutch holder (only in combination with tip assist)	
		Anti-tip tubes swing-away (2)	
KRR090004		Left KRR090005 □ Right (+)	
KRR090006		Pair (+) (not available with tip-assist)	
		7,	
KRR090007		Quickie style anti-tip tubes pluggable (2) (+) (Pair)	
		Anti tip flip up (user operated) (2)	
KRR090033		Left KRR090008 ☐ Right (+) (Please choose only one)	
KRR090010		Transit wheels (2) (+)	
KRR090016		Frame protection, black (textile, with velcro)	
KRR090577		Frame protection, Plastic, glued to the frame, black	
			<u> </u>
		FORKS (compare with your choice on the first page)	
KRR080013	0	Castor fork, Carbotecture ®, 111mm Colour: Black (L)	
		Colour: Silver Orange Red (3) Blue (3)	
KRR080029	0	Castor fork, Aluminium, 123mm Colour: Black	
KBBooooo		Colour: Silver Orange	
KRR080008	П	Castor fork, aluminium, one-arm Colour: Black (only with 4" & 5" Alu soft, Alu solid) Colour: Sliver Orange Red Blue	
KRR080007		Castor fork, aluminium, suspension (Frog Legs ®), black (+)	
7.4.4.4000007		Cade fork, diaminant, cacpondent (110g 20g0 0), black (1)	
		CASTORS (compare with your choice on the first	page)
KRR080301	0	Castor, solid (hard), rim: composite, black (3"-6", not with One arm fork)	
KRR080302	0	Castor, solid (soft), rim: composite, black (only 5" and 6", not with One arm fork)	
KRR080303		Castor, solid (hard), light-up, rim: composite, black (not available in 6", not with One arm fork)	
KRR080304		Castor, solid (hard), rim: aluminium Colour, rim: Silver (only 4" and 5")	
KRR080306		Castor, solid (soft), rim: aluminium Colour, rim: Silver (only 4" and 5")	
		Colour, rim: ☐ Black ☐ Orange ☐ Red ☐ Blue	
		DEAD WHEELS	tuno)
		REAR WHEELS (compare with your choice on the first page; choose only one wheel NOTE: for users with higher user weight than 100 kg, we recommend cross spoked wheels	type)
KRR070000	0	Without rear wheel	
KRR070000 KRR070001	0	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes (+)	
KRR070002	0	Rear wheel, Design, 36 spokes (silver), straight spokes ((nly 24")	
KRR070003	0	Rear wheel, Lightweight, 24 spokes (silver), straight spokes (Std.: black hub/rim; colour options last p	page)
KRR070004		Rear wheel, Proton, 24 spokes (black), straight spokes (L) r weight limit 100kg)	
KRR070005		Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver	
KRR070007		Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub black	
(1):	Not	possible with carbon axle-stems and axle: chair will be built with aluminium axle-stems and axle-tube	For more information
(2):		with KRR07049: camber adapter 0°, superlight.	www.SunriseMedical.eu
(3):		and blue only in combination with soft castor	© Sunrise Medical 21 & 22
	(L)=	lightest option (+)= adds 500 g or more to the wheelchair-weight	

FRAME ACCESSORIES



KRR090032 Suspension



KRR090001 Tip assist



KRR090004 Anti-tip tubes SA



KRR090007 Quickie style anti tip tubes pluggable



KRR090033 Anti tip flip up (user operated)



KRR090010 Transit Wheels



KRR090010 Frame protection,Plastic, glued to the frame, black

FORKS



KRR080013 & KRR080306 Castor fork, Carbotecture & Castor solid (soft) aluminum



KRR080008 Castor fork, aluminum, one-arm, red

REAR WHEELS & HANDRIMS



KRR070001 Universal wheel



KRR070002 Design wheel



KRR070003 Lightweight wheel



KRR070004 Proton wheel



KRR070212 Handrims, Ellipse 3R



KRR070208 Handrims Maxgrepp

TYRES



KRR070101 Schwalbe Right Run



KRR070107 Schwalbe One



KRR070104 Schwalbe Marathon Plus



KRR070102 Tyre, puncture proof

WHEEL LOCKS



KRR060001 Standard brake



KRR060002 Knee lever brake



KRR060023 Compact brake EVO




KRR060024 Lightweight compact brake EVO

Krypton	R		
	•	= Standard O = No cost option \square = Cost option Pr	ice in €
		TYRES (choose only one tyre	type)
KRR070101	0	Tyre, pneumatic, fine profile, Schwalbe Right Run (not in 25")	
KRR070107		Tyre, pneumatic, slick, Schwalbe One (L)	
KRR070102	0	Tyre, puncture proof, tread (+)	
KRR070104		Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution (+)	
KRR090163		Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube (not with solid	tyres)
		HANDRIMS	
KRR070201	0	Handrims, aluminium, powder-coated, clear (L)	
KRR070202	0	Handrims, aluminium, anodised, silver (L)	
KRR070203		Handrims, aluminium, powder-coated, coloured (L) (select colour on last page)	
KRR070204	0	Handrims, stainless steel, silver (+)	
KRR070205		Handrims, titanium, untreated	
KRR070206		Handrims, titanium, powder-coated, clear	
KRR070207		Handrims, Supergrip ® (+)	
KRR070208		Handrims, Maxgrepp ® (+) (only 24")	
KRR070212		Handrims, Ellipse 3R, black (+) (ergonomic handrims with rubber-str	rip)
KRR070250		Handrim cover, rubber (not with 25")	
		Hand rim mounting	
KRR070307		Handrim mount, super narrow (7 mm) (only with stainless steel)	
KRR070317	•	Handrim mount, narrow (17 mm)	
KRR070327	0	Handrim mount, wide (27 mm)	
		MULTI CAMPED (change grity one option)	
KRR070049		WHEEL CAMBER (choose only one option) 0° superlight, no width-adjustment (1) (L)	
KRR070600	0	0° width-adjustment +/-15mm	
KRR070603	0	3° width-adjustment +/-15mm	
		o waar adjustment ii fomiii	
		SPOKE GUARDS (only available in 24") (+)	
KRR070831		Transparent	
KRR07	70813	Galactica	
KRR070823		Princess KRR070818 □ Jungle Fever	
KRR070816		Happy Lions KRR070811 □ Feather	
KRR070819		Monkeys KRR070821 □ Pirates	
KRR070820		Oriental KRR070806 Chrome wheel	
KRR070829 KRR070817		Superhero KRR070814 □ Graffiti Hot wheels KRR070802 □ Basketball KRR07080 □ Basketball 2	
KRR070822		Pool KRR070812 Flying water	
KRR070827		Silk KRR070828 Space	
KRR070830		Tiger KRR070833 □ Zebra	
		WHEEL LOCKS (Choose only one option)	
KRR060001		Standard wheel lock	
KRR060002		Knee lever wheel lock	
KRR060010		Extension for knee lever wheel lock, slip-on, swing-away (only with KRR06)	0002)
KRR060023		Compact wheel lock EVO	
KRR060024		Lightweight compact brake EVO (2) (L)	
KRR060006		One arm wheel lock (left hand actuation) (+)	
KRR060007		One arm wheel lock (right hand actuation) (+)	
		ACCESSORIES	
KRR090000	•	Tool kit	
KRR090024		High pressure pump	
KRR090025		Air pump	
KRR090026		Mobile phone pocket	
KRR090030		Back pack	
			<u> </u>
(1):	Can	ber adapter glued into the axle-tube. No adjustment of the track-width possible.	For more information
(2):	-	stment of COG ± 2 cm on a chair in the field.	www.SunriseMedical.eu
	(L)=	lightest option (+)= adds 500 g or more to the wheelchair-weight	© Sunrise Medical 21 & 22

Krypton	R						
	•		cost option	□ = Cost option			
				HUB/ RIM OF THE STANDARD WH		Hub/rim (4)	
		Handrim (2) H	1ub/IIIfi (1)		Handrim (2)	Hub/rim (1)	
02 Brilliant Silv	er			27 Graphitemetallic			
03 Silver Matt				29 Black			
06 Sun Yellow				31 Matt Black			
14 Blue				32 Signal White			
15 Brilliant Mar	rine			33 Pink			
16 Matt Marine)			34 Dark Chocolate Metallic			
19 Brilliant Pur	ple			35 Skyblue			
21 Red Ruby				36 Cyan-Metallic			
22 Brilliant Red	t			37 Marsred (Orange)			
23 Matt Red				40 Kiwi			
24 Tutti Frutti C	Orar	nge 🗆		41 Flamingo-Glitter			
				42 Matt-Anthracite			
				46 Mamba Green (new)			
(1) Only with ur	niver	sal cross spoked rear wheels		(2) only with KRR070203, handrin	n aluminium, po	wder coated	
		FRAME-COLOURS					
KRR100139	0	Natural Carbon, no coating					
KRR100137		Transparent, gloss					ļ
KRR100138		Transparent, matt					
KRR100122		Black (gloss)					
KRR100123		Matte black					
		FRAME-GRAPHICS "F	lying Pixel"	(scratch resistant, non-removeable, include	s the "Kryptor	n" sticker)	
KRR100590		Frame graphics Flying Pixe		, , , , , , , , , , , , , , , , , , , ,	71		
KRR100591		Frame graphics Flying Pixe	el, red				
KRR100592		Frame graphics Flying Pixe	el, blue				
		COLOUR STICKERS "	SOPUR"/ "Q	QUICKIE" and "Krypton" (3)			
KRR100550	0	black		<i>7</i> 1			
KRR100551	0	silver					*
KRR100552	0	orange					
KRR100553	0	red					,
KRR100554	0	blue					
		COLOURS FOR LIGHT	WEIGHT WI	HEELS 24" (4)			
KRR100336		Rear wheel, lightweight and					
KRR100337		Rear wheel, lightweight and	odized, blue (4)				
KRR100338		Rear wheel, lightweight and	odized, red (4)				
		ole when frame graphics "Flying s for hub and rim, black spokes.					
Camananta	£	anasial buildi					
Comments	TOI	special build:					
					bı	uilt <mark>4</mark> m	
						SOM AUTHORISTICS	D.
Krypton R	₹ -	Base Price					
Option price							
Packaging a		transport					-
Exclusive VA							
VAT.							
Total Pric	е						
Sunrise Medic	cal (Customer Service				For more informatio	
		chring 2-4, D-69254 Malsch	/Heidelberg	E-mail: export@sunrismedical.d		ww.SunriseMedical	
Tel. +497253				Website: sunrisemedical.eu	© Su	nrise Medical 2021	1 & 22
Fax. +497253	980	-222					

FRAME COLOURS



- Transparent, gloss Transparent, matt Black, gloss Black, matt

FRAME GRAPHIC "FLYING PIXELS"



From left to right:

- Frame graphics Flying Pixel, red Frame graphics Flying Pixel, silver Frame graphics Flying Pixel, blue

24 " LIGHT WEIGHT WHEEL COLOURS



KRR100336 Lightweight wheel, anodized orange



KRR100337 Lightweight wheel, anodized blue



KRR100338 Lightweight wheel, anodized red