



Deliver to:  
 Name:  
 Street:  
 Town: Postal code:  
 E-Mail:  
 Tel: Fax:



Invoice to:  
 Name:  
 Street: Postal code:  
 Town: Postal code:  
 E-Mail:  
 Tel: Fax:

Order date:  
 Purchase order:  
 Number required:

Marked for:  
 Sunrise order no:

- Order    Calculation    Built 4 me

- = Standard    = Cost option

- = No cost option

**Customer service fax number: +497253 980-222**

**Max. user weight: 140 kg <sup>(1)</sup>**

**Reimbursement code:**

**USER WEIGHT**  kg

**Art. no.: 05300001 IC**

Mandatory information to ensure correct configuration

**Base Price**

**CLIENT MEASUREMENTS** (all measurements without cushion)

**SEAT WIDTH SELECTION** - measured from outside frame tube to outside frame tube (mm)

300	320	340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Includes option NSW016602:  
 Userweight 140 kg,

**SEAT DEPTH SELECTION** - in mm

340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

(1) Max user weight 140kg only with option NSW010602 or with SW 440-500 mm;  
 Standard maximum user weight: 120kg.

● = Standard    ○ = No cost option    □ = Cost option

## FRONT SEAT HEIGHT SELECTION: in dependance of front castors, forks and frame versions

		Front seat height in mm																								
Frame		Fork (length)	350	360	370	380	390	400	410	420	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570	
Front castor	3"	111 mm									○	○	○	○	○											
		123 mm* (1)											○	○	○	○										
		138 mm (1)												○	○	○	○									
		174 mm (1)															○	○	○	○	○					
		One Arm (1)										○	○	○	○	○										
	Hemi	111 mm* (1)	□	□	□																					
		138 mm* (1)		□	□	□	□																			
		174 mm* (1)					□	□	□	□	□															
	4"	111 mm										○	○	○	○	○										
		123 mm* (1)												○	○	○	○									
		138 mm (1)													○	○	○	○	○							
		174 mm (1)															○	○	○	○	○					
		One Arm (1)											○	○	○	○	○									
	Hemi	111 mm* (1)	□	□	□	□																				
		138 mm* (1)			□	□	□	□	□																	
		174 mm* (1)					□	□	□	□	□	□														
	5"	111 mm													○	○	○	○								
		123 mm* (1)														○	○	○	○							
138 mm (1)																○	○	○	○	○						
174 mm (1)																	○	○	○	○	○	○	○			
One Arm (1)																○	○	○	○							
Hemi	111 mm* (1)			□	□	□	□																			
	138 mm* (1)				□	□	□	□	□																	
	174 mm* (1)					□	□	□	□	□	□															
6"	111 mm															○	○									
	123 mm* (1)																○	○	○							
	138 mm (1)																○	○	○	○	○					
	174 mm (1)																	○	○	○	○	○	○	○		
	One Arm (1)																○	○								
Hemi	111 mm* (1)				□	□																				
	138 mm* (1)					□	□	□	□	□																
	174 mm* (1)						□	□	□	□	□	□														
7"	123 mm* (1)																		○	○						
	138 mm (1)																			○	○	○				
	174 mm (1)																				○	○	○	○	○	
	Hemi 174 mm* (1)											□	□	□	□	□	□									

\* only in black colour available

## REAR SEAT HEIGHT SELECTION: in dependance of rear wheel size

		Rear seat height in mm																				
Rear wheels	330	340	350	360	370	380	390	400	410	420	430	440	450	460	470	480	490	500	510	520	530	
22"	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○						
24"				○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○		
25"					○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
26"						○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

(1) Compatible with power add on devices.

For more information  
visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 21 & 22

### Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

#### Seat width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether)
- b) On the wheelchair: From outside frame tube to outside frame tube.

#### Seat depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee. (approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

#### Front seat height:

- a) For the user: When sitting with 90° kneeangle: From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground. (min. 30 mm longer than user measurement) Remember to include the cushion height.

#### Rear seat height:

- b) On the wheelchair: From back of seat tube to the ground.  
Note: Max. seat height difference rear to front is 110 mm

#### Back height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.  
Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably

#### Back angle:

Important: The back height depends on the degree of disability.  
0° is the angle 90° to the ground. Negative angles are <90° and positive angles are >90°.

#### Technical Data Overview:

<b>Max. overall height:</b>	1020 mm
<b>Max. overall length:</b>	1150 mm
<b>Max. pivot width:</b>	700 mm
<b>Max. safe slope:</b>	7° (for use of parking brakes)
<b>Wheelchair weight:</b>	from 11,5 kg
<b>Dismantled weight:</b>	from 9,1 kg
<b>Mass of heaviest part:</b>	2,1 kg
<b>Max. overall width:</b>	<b>Camber</b>
	For 0° camber: SW + 160 mm
	For 2° camber: SW + 220 mm
	For 4° camber: SW + 270 mm

**Intended use and environment** - for people with limited mobility; self-propelling or being pushed by an assistant, indoor and outdoor use.

Product is **crash tested according to ISO 7176-19** and therefore approved for transportation in a vehicle in combination with option "Tie down brackets for transportation". For transport without a user fix the chair safe in the vehicle (e.g. car or plane)

All changes in the configuration might change the seat heights.  
Manufacturing tolerances of +/- 10 mm

# Neon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## SEAT

### Seat sling

NSW020002	●	Seat sling, nylon, tension adjustable, 1 utility bag	
NSW020003	□	Seat sling, nylon, tension adjustable, 2 utility bags	(only from SW 380 mm)
NSW020004	□	Seat sling, nylon, tension adjustable, 1 utility bag, 1 catheter bag	(only from SW 340 mm)
NSW020000	○	Without seat sling	

### Seat cushion, cover with zipper, black

Foam (medium-soft)	Thickness:	NSW020101	□	30 mm	NSW020102	□	50 mm
Latex (very soft)	Thickness:	NSW020103	□	50 mm			

### Seat accessories

NSW020050	□	Additional padding for seat upholstery, removable	
NSW090018	□	Lap belt with buckle, metal	

### JAY Cushions - for complete product range see Jay order forms

JAY020002	□	JAY Basic	
JAY020003	□	JAY Soft Combi P	
JAY020001	□	JAY Easy Visco, Curved Base	
JAY020004	□	JAY Lite	
JAY020007	□	JAY Easy Fluid, Curved Base	

## BACKREST

### Back height range (mm)

Choose back height factory setting.

Height adjustments of in total 50 mm possible.

250 - 300	300 - 350	350 - 400	400 - 450	450 - 475
○ 250	○ 300	○ 350	○ 400	○ 450
○ 275	○ 325	○ 375	○ 425	○ 475
○ 300	○ 350	○ 400	○ 450	

NSW030010	○	Backrest, fixed (90°)	
NSW030031	□	Backrest, fold-down (to the front)	(1)
NSW030012	□	Backrest, angle-adjustable	(75° to 103°)
NSW030031	□	Backrest, fold-down (to the front)	(1)

### Pushhandles

NSW030200	○	Without push handles	
NSW030201	○	Push handles long	
NSW030203	□	Push handles, fold-down	(available from back height 325 mm)
NSW030205	□	Push handles, height-adjustable for Jay backs	(available from back height 300 mm)
NSW030204	□	Push handles, height-adjustable	(available from back height 300 mm)

### Backrest upholstery

NSW030316	○	EXO Backrest sling, padded, adjustable, 1 accessory bag, black (EXO Evo version)	
		Colour, binding tape: ● Black    ○ Orange    ○ Silver    ○ Red    ○ Blue	
NSW030317	□	EXO PRO Backrest, padded, adj., breathable spacer fabric, 1 accessory bag, black (EXO Evo version)	
		Colour, binding tape: ● Black    ○ Orange    ○ Silver    ○ Red    ○ Blue	
NSW030300	○	Without backrest upholstery	


### Accessories for backrest

NSW030402	□	Stabilising bar, automatic folding	(only available from SW 360- 480 mm)
NSW030401	□	Stabilising bar, swing-away	
NSW030412	□	Headrest, medium, height, depth and angle-adjustable	
NSW030411	□	Headrest, tall, height, depth and angle-adjustable, 200x110 mm	

(1)

Not with height adjustable push handles, JAY backs, stabilising bar or headrest. Please note folding-Matrix on next page. If push handles long chosen only from back height 300 mm available; if push handles fold-down chosen only from back height 350 mm available

For more information  
visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 21 & 22



NSW020002  
Seat sling with 1 utility bag



NSW020003  
Seat sling with 2 utility bags



NSW090018  
Lap belt with buckle, metal



NSW030012  
Backrest, angle adjustable



NSW030201  
Push handles long



NSW030203  
Push handles, fold-down



NSW030204  
Push handles, height-adjustable



NSW030316  
EXO backrest sling



NSW030317  
EXO PRO backrest sling



NSW030412  
Headrest, medium

**Folding-Matrix for NSW030031 Backrest, fold-down to the front.**



NSW030031  
Backrest, fold-down

		Back height						
		325	350	375	400	425	450	475
Seat width	320	✓	✓	✓	✓	✓	✓	✓
	340	✓	✓	✓	✓	✓	✓	✓
	360	✓	✓	✓	✓	✓	✓	✓
	380		✓	✓	✓	✓	✓	✓
	400			✓	✓	✓	✓	✓
	420				✓	✓	✓	✓
	440				✓	✓	✓	✓
	460					✓	✓	✓
	480						✓	✓
	500							✓



=folds 90° to the front on folded chair



=can't be folded full 90° to the front on folded chair in combination with fold-down pushhandles



=can't be folded full 90° to the front on folded chair



=will fold less than 45° to the front on folded chair, folds 90° on unfolded chair

# Neon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## JAY J3 BACKS

### Jay J3 Back - Shallow Contour

Mid Thoracic:	<input type="checkbox"/> JAY030010	Medium - 340 mm	<input type="checkbox"/> JAY030033	Tall - 380 mm
Upper Thoracic:	<input type="checkbox"/> JAY030012	Short - 420 mm	<input type="checkbox"/> JAY030019	Medium - 460 mm
Shoulder Height:	<input type="checkbox"/> JAY030015	Short - 530 mm	<input type="checkbox"/> JAY030022	Medium - 570 mm

### Jay J3 Back - Mid Contour

Mid Thoracic:	<input type="checkbox"/> JAY030011	Medium - 340 mm	<input type="checkbox"/> JAY030018	Tall - 380 mm
Upper Thoracic:	<input type="checkbox"/> JAY030013	Short - 420 mm	<input type="checkbox"/> JAY030020	Medium - 460 mm
Shoulder Height:	<input type="checkbox"/> JAY030016	Short - 530 mm	<input type="checkbox"/> JAY030023	Medium - 570 mm

### JAY 3 Back - Mid Deep Contour (150 mm Contour Depth)

Mid Thoracic:	<input type="checkbox"/> JAY030044	Medium - 340 mm	<input type="checkbox"/> JAY030041	Tall - 380 mm
Upper Thoracic:	<input type="checkbox"/> JAY030045	Short - 420 mm	<input type="checkbox"/> JAY030042	Medium - 460 mm
Shoulder Height:	<input type="checkbox"/> JAY030043	Short - 530 mm	<input type="checkbox"/> JAY030047	Medium - 570 mm

### Jay J3 Carbon Back - Shallow Contour

Lower Thoracic:	<input type="checkbox"/> JAY060001	Short - 170 mm	<input type="checkbox"/> JAY060002	Tall - 240 mm
Mid Thoracic:	<input type="checkbox"/> JAY060003	Short - 300 mm	<input type="checkbox"/> JAY060004	Tall - 380 mm

### Accessories for JAY J3 backs

WHT030415	<input type="checkbox"/> Whitmyer Headrest, COBRA XTRA, flip back	<i>(only with JAY J3 backs)</i>
WHT030416	<input type="checkbox"/> Whitmyer Headrest, AXYS, Standard	<i>(only with JAY J3 backs)</i>
JAY090058	<input type="checkbox"/> JAY Accessory bag	<i>(only with Jay J3 Carbon backs)</i>

## ARMREST AND SIDEGUARDS

### Armrest

NSW040005	<input type="checkbox"/> Desk sideguard, flip-up, removable	
NSW040001	<input type="checkbox"/> Desk sideguard, armpad (250 mm), fixed, flip-up, removable	
NSW040002	<input type="checkbox"/> Desk sideguard, armpad (250 mm), height-adjustable, flip-up, removable	
NSW040003	<input type="checkbox"/> Desk sideguard, armpad (310 mm), fixed, flip-up, removable	
NSW040004	<input type="checkbox"/> Desk sideguard, armpad (310 mm), height-adjustable, flip-up, removable	
NSW040050	<input type="checkbox"/> Single-post sideguard, armpad (250 mm), height-adjustable, removable <b>(1)</b>	
NSW040051	<input type="checkbox"/> Single-post sideguard, armpad (360 mm), height-adjustable, removable <b>(1)</b>	
NSW040052	<input type="checkbox"/> Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable <b>(1)</b>	
NSW040117	<input type="checkbox"/> Tubular armrest, Aluminium, swing away (sideways), padded, removable	<i>(only with backrest, angle adj.)</i>

### Accessories and options for desk sideguards

NSW040020	<input type="checkbox"/> Release assistance for tetraplegics	<i>(only with desk sideguard)</i>
NSW090012	<input type="checkbox"/> Tray, swing-away (to the side), SW 340 - 500 mm, left-mount	<i>(only with desk sideguard)</i>
NSW090013	<input type="checkbox"/> Tray, swing-away (to the side), SW 340 - 500 mm, right-mount	<i>(only with desk sideguard)</i>
NSW090011	<input type="checkbox"/> Tray, push-fit, SW 380 - 480 mm	<i>(only with desk sideguard)</i>

### Sideguards

NSW040107	<input type="checkbox"/> Sideguard, composite, removable, 210 mm high (310 mm long), black	
NSW040102	<input type="checkbox"/> Sideguard, aluminium, cold weather protection (in frame colour)	
NSW040103	<input type="checkbox"/> Sideguard, aluminium, fender, cold weather protection (in frame colour)	
NSW040104	<input type="checkbox"/> Sideguard, carbon fibre, fender, cold weather protection (black)	

## HANGER

NSW050001	<input type="checkbox"/> Hanger angle 80°	
NSW050002	<input type="checkbox"/> Hanger angle 70° <b>(2)</b>	

### Hanger, aluminium, swing-away

NSW050006	<input type="checkbox"/> Hanger, aluminium, swing-away, LH	
NSW050007	<input type="checkbox"/> Hanger, aluminium, swing-away, RH	
NSW050008	<input type="checkbox"/> Hanger, aluminium, swing-away, pair	

### Hanger, aluminium, flip-up (0-90°), swing-away

NSW050009	<input type="checkbox"/> Hanger, aluminium, flip-up (0° - 90°), swing-away, LH	
NSW050010	<input type="checkbox"/> Hanger, aluminium, flip-up (0° - 90°), swing-away, RH	
NSW050011	<input type="checkbox"/> Hanger, aluminium, flip-up (0° - 90°), swing-away, pair	<i>(from SW 380 mm)</i>

- (1) Only with COG Adapter Standard  
 (2) In combination with NSW010024 "Frame, Hemi" only 70° hangers available

For more information  
 visit our webpage

[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

© Sunrise Medical 2021 & 22

JAY J3 Back  
Shallow Contour



JAY J3 Back  
Mid Contour



JAY J3 Back  
Mid Deep Contour



Jay 3 Carbon Back  
Shallow Contour



NSW04000x  
Desk sideguard, armpad, flip-up, removable



NSW040050  
Single-post sideguard




NSW040052  
Single-post sideguard, armpad (330 mm)



NSW040117  
Tubular armrest, swing away, removable, padded



NSW040107  
Sideguard, composite, removable (black)



NSW040103  
Aluminium, cold weather protection, fender



NSW040104  
Carbon fibre, cold weather protection, fender



NSW050011  
Hanger, aluminum, flip-up (0-90°), Swing-away

# Neon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## FOOTBOARDS

NSW050104	<input type="checkbox"/>	Footboard, divided, composit, flip-up, heel loop, black	(from SW 380 mm)	
NSW050020	<input type="checkbox"/>	Footboard, divided, composit, angle-adjustable, flip-up, heel loop, black	(from SW 380 mm)	
NSW050022	<input type="checkbox"/>	Footboard, divided, aluminum, angle-adjustable, flip-up, heel loop, frame colour		
NSW050051	<input type="checkbox"/>	Footboard, divided, aluminum, angle-adjustable, flip-up, high-mount, heel loop, frame colour	(1)	
NSW050023	<input type="checkbox"/>	Footboard, platform, aluminum, angle-adjustable, flip-up, heel loop, frame colour		
NSW050126	<input type="checkbox"/>	Swing-away mechanism for platform footboard		
NSW050059	<input type="checkbox"/>	Footboard, platform, Performance, aluminum, angle-adjustable, flip-up, calf strap	(from SW 320 mm)	
NSW050131	<input type="checkbox"/>	Footboard, platform, lightweight, plastic, angle & depth-adjustable flip up, calf strap, black		
NSW050132	<input type="checkbox"/>	Footboard, platform, lightweight, carbon fibre, angle & depth-adjustable, flip up, calf strap, black	(1)	

## Accessories for footboards

NSW050127	<input type="checkbox"/>	Adjustable side protection, to create a side-stop on the footboard	(2)	
-----------	--------------------------	--	-----	--

## LOWER LEG LENGTH (LLL)

(LLL has to be 25 mm less than front seat height; measurements in mm)

- 220    230    240    250    260    270    280    290    300    310  
 320    330    340    350    360    370    380    390    400    410  
 420    430    440    450    460    470    480    490    500    510

LLL below approximately 320 mm require high-mount footboards- depends on configuration, please check configurator for details

## FRAME VERSIONS

NSW010015	<input type="radio"/>	Frame SA		
NSW010602	<input type="checkbox"/>	Frame reinforcement for max. user weight 140 kg		
NSW010024	<input type="radio"/>	Frame, Hemi for lower Seat heights	(compare with your choice on the first page)	

## CENTRE OF GRAVITY (COG)

(20 mm= very passive; 120 mm= very active)

NSW070703	Adapter standard	<input type="radio"/> 20	<input type="radio"/> 30	<input type="radio"/> 40	<input type="radio"/> 50	<input type="radio"/> 60	<input type="radio"/> 70	<input type="radio"/> 80	<input type="radio"/> 90
		<input type="radio"/> 100	<input type="radio"/> 110	<input type="radio"/> 120					
NSW070701	Adapter lightweight	<input type="radio"/> 50	<input type="radio"/> 60	<input type="radio"/> 70	<input type="radio"/> 80	<input type="radio"/> 90	<input type="radio"/> 100	<input type="radio"/> 110	<input type="radio"/> 120

## FRAME ACCESSORIES

(choose either tip assist or anti- tip tube per side)

NSW090016	<input type="checkbox"/>	Frame protection, black		
NSW090100	<input type="checkbox"/>	Tie down brackets for transportation		
NSW090010	<input type="checkbox"/>	Transit wheels	(not in combination with anti-tip tubes or tip assist)	

## Anti-tip tubes, swing-away

NSW090004	<input type="checkbox"/>	Left	NSW090005	<input type="checkbox"/>	Right
NSW090006	<input type="checkbox"/>	Pair			(no additional tip assist possible)

## Tip assist

NSW090001	<input type="checkbox"/>	Left	NSW090002	<input type="checkbox"/>	Right
NSW090003	<input type="checkbox"/>	Tip assist, pair			(no additional anti-tip tube possible)
NSW090009	<input type="checkbox"/>	Crutch holder, loop			(only in combination with tip assist)

## EMPULSE F55 BUNDLING

\* Comming soon

EMO000101	<input type="radio"/>	With F55 order, not mounted to the wheelchair	(Please use separate OF to order the F55.)	
EMO000102	<input type="radio"/>	With F55 order, mounted to the wheelchair	(Please use separate OF to order the F55.)	

## FORKS

(compare with your choice on the first page)

NSW080013	<input checked="" type="radio"/>	Castor fork, Carbotecture®, 111mm	Colour: Black	
		Colour:	<input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red (3) <input type="checkbox"/> Blue (3)	
NSW080029	<input type="radio"/>	Castor fork, Aluminium	Colour: Black	(only for 123mm, 138mm and 174mm forks)
		Colour:	<input type="checkbox"/> Silver <input type="checkbox"/> Orange	
NSW080008	<input type="checkbox"/>	Castor fork, aluminium, one-arm	Colour: Black	
		Colour:	<input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue	

## STANDARD CASTORS

(compare with your choice on the first page)

NSW080301	<input type="radio"/>	Castor, solid (hard), rim: composite, black		
NSW080302	<input type="radio"/>	Castor, solid (soft), rim: composite, black		(only 5" & 6")
NSW080303	<input type="checkbox"/>	Castor, solid (hard), light-up, rim: composite, black		(only in 4" and 5")
NSW080304	<input type="checkbox"/>	Castor, solid (hard), rim: aluminium	Colour, rim: Silver	(only 4" & 5")
NSW080305	<input type="radio"/>	Castor, pneumatic, rim: composite, black		(only 7")
NSW080306	<input type="checkbox"/>	Castor, solid (soft), rim: aluminium	Colour, rim: Silver	(only 4" & 5")
		Colour, rim:	<input type="checkbox"/> Black <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue	

- (1)** From SW 340 mm available  
**(2)** Only for footboard, platform lightweight or Performance  
**(3)** Only available in combination with castor, soft (NSW080302 or NSW080306)

For more information  
 visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
 © Sunrise Medical 21 & 22





NSW050022  
Footboard divided,  
aluminium, frame colour



NSW050059  
Footboard platform,  
Performance



NSW050132  
Footboard platform,  
lightweight, carbon fibre



NSW050127  
Adjustable side  
protection



NSW090010  
Transit wheels



NSW090006  
Anti-tip tubes, swing away



NSW090001  
Tip assist



NSW080013 & NSW080306  
Castor fork, Carbotecture & Castor solid (soft) aluminum



NSW080008  
Castor fork, aluminum, one-arm, red

**F55 BUNDLING - NOGO MATRIX FOR 24" BACKWHEELS**

COG	Seat Depth						
	380	400	420	440	460	480	500
40	X	X	X	✓	✓	✓	✓
60	X	X	X	X	✓	✓	✓
70	X	X	X	X	✓	✓	✓
80	X	X	X	X	X	✓	✓
90	X	X	X	X	X	✓	✓
100	X	X	X	X	X	X	✓
110	X	X	X	X	X	X	✓
120	X	X	X	X	X	X	X
130	X	X	X	X	X	X	X

# Neon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## REAR WHEELS

(compare with your choice on the first page)

NSW070000	<input type="radio"/>	Without rear wheel	
NSW070001	<input type="radio"/>	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes	(not in 25")
NSW070002	<input type="radio"/>	Rear wheel, Design, 36 spokes (silver), straight spokes	(only 24" and 26")
NSW070003	<input type="radio"/>	Rear wheel, Lightweight, 24 spokes (silver), straight spokes	(not in 26"; see colour options on last page)
NSW070004	<input type="checkbox"/>	Rear wheel, Proton, 24 spokes (black), straight spokes	(not in 26"; user weight limit 100kg)
NSW070005	<input type="checkbox"/>	Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver	(not in 22")
NSW070006	<input type="checkbox"/>	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub: silver	(only 24" & 25")
NSW070007	<input type="checkbox"/>	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub black	(only 24" & 25")
NSW070008	<input type="checkbox"/>	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes	(1) (only 22" & 24")
NSW070009	<input type="checkbox"/>	Rear wheel, Mountainbike, 36 spokes (silver), straight spokes	(only 24")
NSW070010	<input type="checkbox"/>	Mountainbike wheel, replacement kit incl. tyre (extra coarse profile)	(only 24")

## Axle

NSW070150	<input checked="" type="radio"/>	Quick-release axle	
NSW070151	<input type="checkbox"/>	Quick-release axle, designed for quads	(2)

## TYRES

NSW070101	<input type="radio"/>	Tyre, pneumatic, fine profile, Schwalbe Right Run	(not in 25")
NSW070102	<input type="radio"/>	Tyre, puncture proof, tread	(not in 26")
NSW070106	<input type="checkbox"/>	Tyre, pneumatic, tread, Schwalbe Downtown Two Grip	(only 24")
NSW070107	<input type="checkbox"/>	Tyre, pneumatic, slick, Schwalbe One	(only 24" and 25")
NSW070104	<input type="checkbox"/>	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	

NSW090163	<input type="checkbox"/>	Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube	(not with solid tyres)
-----------	--------------------------	--	------------------------

## HANDRIMS

NSW070204	<input type="radio"/>	Handrims, stainless steel, silver	
NSW070202	<input type="radio"/>	Handrims, aluminium, anodised, silver	
NSW070201	<input type="radio"/>	Handrims, aluminium, powder-coated, clear	
NSW070203	<input type="checkbox"/>	Handrims, aluminium, powder-coated, coloured	(select colour on last page)
NSW070205	<input type="checkbox"/>	Handrims, titanium, untreated	(not in 22")
NSW070206	<input type="checkbox"/>	Handrims, titanium, powder-coated, clear	(not in 22")
NSW070207	<input type="checkbox"/>	Handrims, Supergrip ®	
NSW070208	<input type="checkbox"/>	Handrims, Maxgripp ®	(only 22" & 24")
NSW070212	<input type="checkbox"/>	Handrims, Ellipse 3R, black	(ergonomic handrims with rubber-strip, 24" and 25" only)

## Hand rim mounting

NSW070307	<input type="checkbox"/>	Handrim mount, super narrow (7 mm)	(only with stainless steel)
NSW070317	<input checked="" type="radio"/>	Handrim mount, narrow (17 mm)	
NSW070327	<input type="radio"/>	Handrim mount, wide (27 mm)	

## Handrim accessories

NSW070250	<input type="checkbox"/>	Handrim cover, rubber	(not with 25" and 26")
-----------	--------------------------	-----------------------	------------------------

## WHEEL CAMBER

NSW070600	<input type="radio"/>	0°	
NSW070602	<input type="radio"/>	2°	(3)
NSW070604	<input type="radio"/>	4°	(3)

## SPOKE GUARDS

(coloured: only 22 & 24"; transparent: only 22", 24" & 26")

NSW070831	<input type="checkbox"/>	Transparent						
		NSW070813	<input type="checkbox"/>	Glaktica				
NSW070823	<input type="checkbox"/>	Princess	NSW070818	<input type="checkbox"/>	Jungle Fever			
NSW070816	<input type="checkbox"/>	Happy Lions	NSW070811	<input type="checkbox"/>	Feather			
NSW070819	<input type="checkbox"/>	Monkeys	NSW070821	<input type="checkbox"/>	Pirates			
NSW070820	<input type="checkbox"/>	Oriental	NSW070811	<input type="checkbox"/>	Feather			
NSW070829	<input type="checkbox"/>	Superhero	NSW070814	<input type="checkbox"/>	Graffiti			
NSW070817	<input type="checkbox"/>	Hot Wheels	NSW070802	<input type="checkbox"/>	Basketball	NSW070803	<input type="checkbox"/>	Basketball 2
NSW070822	<input type="checkbox"/>	Pool	NSW070812	<input type="checkbox"/>	Flying Water			
NSW070827	<input type="checkbox"/>	Silk	NSW070828	<input type="checkbox"/>	Space			
NSW070830	<input type="checkbox"/>	Tiger	NSW070833	<input type="checkbox"/>	Zebra			

- (1) Not with push handles fold down & not without push handles  
only with adapter, standard and only with 0° and 4° wheel camber
- (2) Not with rear wheels Spinergy / mountainbike/ drum brakes
- (3) Not with COG Adapter, lightweight

For more information  
visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 21 & 22



NSW070001  
Universal wheel



NSW070002  
Design wheel



NSW070003  
Lightweight wheel



NSW070004  
Proton wheel



NSW070010  
Mountainbike wheel



NSW070101  
Schwalbe Right Run



NSW070107  
Schwalbe One



NSW070104  
Schwalbe Marathon Plus  
Evolution



NSW070102  
Tyre, puncture proof



NSW070208  
Handrims Maxgrepp



NSW070212  
Handrims Ellipse 3R

# Neon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## WHEEL LOCKS

(choose only one option)

NSW060001	● Standard brake		
NSW060002	□ Knee-lever brake		
NSW060010	□ Brake lever extension, knee-lever brake, fold-down		
NSW060015	□ Push to lock brake, one-arm, left actuation		
NSW060016	□ Push to lock brake, one-arm, right actuation		
NSW060003	□ Compact brake		

## ACCESSORIES

NSW090000	● Tool kit		
NSW090025	□ Air pump		
NSW090024	□ High pressure pump (0 - 11 bar)		
NSW090030	□ Back pack		
NSW090026	□ Mobile phone pocket		

## COLOURS

(Please consult our colour catalogue to review your colour choice)

	Frame	Handrim	Hub/rim (1)		Frame	Handrim	Hub/rim (1)
02 Brilliant Silver	○	□	□	27 Graphitemetallic	○	□	□
03 Silver Matt	○	□	□		○	□	□
	○	□	□	29 Black	○	□	□
06 Sun Yellow	○	□	□	31 Matt Black	○	□	□
	○	□	□	32 Signal White	○	□	□
	○	□	□	33 Pink	○	□	□
14 Blue	○	□	□	34 Dark Chocolate Metallic	○	□	□
15 Brilliant Marine	○	□	□	35 Skyblue	○	□	□
16 Matt Marine	○	□	□	36 Cyan-Metallic	○	□	□
19 Brilliant Purple	○	□	□	37 Marsred (Orange)	○	□	□
21 Red Ruby	○	□	□	40 Kiwi	○	□	□
22 Brilliant Red	○	□	□	41 Flamingo-Glitter	○	□	□
23 Matt Red	○	□	□	42 Matt Anthracite	○	□	□
24 Tutti Frutti Orange	○	□	□	46 Mamba Green (new)	○	□	□

(1) only with rear wheel, Universal 22" and 24"

## COLOURS FOR REAR WHEEL, LIGHTWEIGHT, 24"

(hub and rim anodized, black spokes)

NSW100336	□ Anodized, orange	(only with Rear wheel, lightweight in 24")	
NSW100337	□ Anodized, blue	(only with Rear wheel, lightweight in 24")	
NSW100338	□ Anodized, red	(only with Rear wheel, lightweight in 24")	

## COLOURS - for stickers ("Neon<sup>2</sup>" on the frame and "Quickie" on sideguard)

NSW100550	○ black		
NSW100551	○ silver		
NSW100552	○ orange		
NSW100553	○ red		
NSW100554	○ blue		

Comments for special build - please apply for a quote; lead time conditions my change.



## Quickie Neon<sup>2</sup> - Basic Price

Option price	.....	€
Packaging and transport	.....	€
Exclusive VAT	.....	€
VAT.	.....	€
<b>Total Price</b>	.....	€

Sunrise Medical Customer service  
 Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg  
 Tel. +497253 980-0  
 Fax. +497253 980-222

E-mail: [export@sunrisemedical.de](mailto:export@sunrisemedical.de)  
 Website: [sunrisemedical.eu](http://sunrisemedical.eu)

For more information  
 visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
 © Sunrise Medical 2021 & 22



NSW060001  
Standard brake



NSW060002  
Knee lever brake



NSW060003  
Compact brake



NSW090000  
Tool kit



NSW090024  
High pressure pump



NSW090025  
Air pump



NSW090030  
Back pack



NSW090026  
Mobile phone pocket



NSW100336  
Lightweight wheel, anodized orange



NSW100337  
Lightweight wheel, anodized blue



NSW100338  
Lightweight wheel, anodized red